

Preparation of genomic DNA from mouse ear punches

1. Obtain one ear punch, place in eppendorf tube.
2. Add 100 μ l ear buffer: 10 mM Tris, pH 8.0 + 50 mM KCl +
0.5% NP40 + 0.5% Tween
3. Heat to 95° C for 15 minutes.
4. Add 2 μ l proteinase K (stock: 10 mg/ml).
5. Incubate at 56° C for 3-4 hours. Vortex every hour or so.
6. Heat to 95° C for 15-20 minutes.
7. (Freeze samples at -20° C.)

PCR using Invitrogen HotWax beads for "hot-start"

1. Prepare pre-mix of all reaction components except DNA and MgCl₂:
 - 5 μ l 5X PCR buffer (Invitrogen; pH 8.5, 9.0, 9.5, or 10.0; see below)
 - 5 μ l dNTP mix (stock: 2 μ M each dNTP)
 - 23-28 μ l H₂O (calculated to have final reaction volume reach 50 μ l)
 - 5 μ l each oligonucleotide (working stock: 10 μ M)
 - 0.5 μ l Taq polymerasemix all of these components in multiple needed (e.g., to test 8 samples make 10X pre-mix by mixing 50 μ l PCR buffer with 50 μ l dNTPs, etc.)
2. Thaw genomic DNA, heat to 95° C 10 min., place on ice.
3. Aliquot 48 μ l pre-mix per geneamp or microamp reaction tube.
4. Add 2 μ l genomic DNA to tubes containing pre-mix.
5. Add one Hotwax bead containing MgCl₂ (1.5, 2.5, or 3.5 mM MgCl₂)
6. Begin PCR cycles:
 - typical settings: 95° C 1 min./55° C 1 min./72° C 2 min. 35X (file #18 on 4th floor thermocycler) or 94° C 45 sec./55° C 1:30/72° C 2 min. 35X (file #109 on 5th floor thermocycler).

The hotwax beads will melt and release the MgCl₂ needed in the reaction only after reaching 95° C; it is unnecessary to add mineral oil.

For novel reactions, optimize by testing a range of rxn buffers (pH 8.5, 9.0, 9.5 or 10.0) versus a range of concentrations of MgCl₂ (1.5, 2.5, or 3.5 mM).